

Anti-ACBD6 antibody (1-282) (STJ118358)

STJ118358

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-ACBD6 (1-282) is suitable for use in Western Blot.
Applications	WB
Host/Source	Rabbit
Reactivity	Human, Mouse, Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	WB 1:500-1:2000
Formulation	PBS containing 0.02% Sodium Azide, 50% Glycerol, pH7.3.
Isotype	IgG
Storage Instruction	Store in a freezer at -20°C and avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	84320
Gene Symbol	ACBD6
Uniprot ID	ACBD6_HUMAN
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-282 of human ACBD6 (NP_115736.1).
Immunogen Region	1-282
Specificity	
Immunogen Sequence	