

Anti-SH2D1B antibody (1-132) (STJ118451)

STJ118451

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-SH2D1B (1-132) is suitable for use in Immunofluorescence.
Applications	IF
Host/Source	Rabbit
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	IF 1:50-1:200
Formulation	PBS containing 0.02% Sodium Azide, 50% Glycerol, pH7.3.
Isotype	IgG
Storage Instruction	Store in a freezer at -20°C and avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	117157
Gene Symbol	SH2D1B
Uniprot ID	SH21B_HUMAN
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-132 of human SH2D1B (NP_444512.2).
Immunogen Region	1-132
Specificity	
Immunogen Sequence	