

Anti-PROZ antibody (220-340) (STJ118523)

STJ118523

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-PROZ (220-340) is suitable for use in Western Blot.
Applications	WB
Host/Source	Rabbit
Reactivity	Mouse

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	WB 1:500-1:2000
Formulation	PBS containing 0.02% Sodium Azide, 50% Glycerol, pH7.3.
Isotype	IgG
Storage Instruction	Store in a freezer at -20°C and avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	8858
Gene Symbol	PROZ
Uniprot ID	PROZ_HUMAN
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 220-340 of human PROZ (NP_003882.1).
Immunogen Region	220-340
Specificity	
Immunogen Sequence	