

Anti-STAP2 antibody (130-220) (STJ118570)

STJ118570

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-STAP2 (130-220) is suitable for use in Western Blot.
Applications	WB
Host/Source	Rabbit
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	WB 1:500-1:2000
Formulation	PBS containing 0.02% Sodium Azide, 50% Glycerol, pH7.3.
Isotype	IgG
Storage Instruction	Store in a freezer at -20°C and avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	55620
Gene Symbol	STAP2
Uniprot ID	STAP2_HUMAN
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 130-220 of human STAP2 (NP_001013863).
Immunogen Region	130-220
Specificity	
Immunogen Sequence	