

## Anti-GPR155 antibody (721-870) (STJ118613)

STJ118613

### GENERAL INFORMATION


|                          |                                                                                       |
|--------------------------|---------------------------------------------------------------------------------------|
| <b>Product Type</b>      | Primary antibodies                                                                    |
| <b>Short Description</b> | Rabbit polyclonal antibody anti-GPR155 (721-870) is suitable for use in Western Blot. |
| <b>Applications</b>      | WB                                                                                    |
| <b>Host/Source</b>       | Rabbit                                                                                |
| <b>Reactivity</b>        | Mouse, Rat                                                                            |

### PRODUCT PROPERTIES

|                            |                                                           |
|----------------------------|-----------------------------------------------------------|
| <b>Clonality</b>           | Polyclonal                                                |
| <b>Clone ID</b>            |                                                           |
| <b>Concentration</b>       |                                                           |
| <b>Conjugation</b>         | Unconjugated                                              |
| <b>Purification</b>        | Affinity purification                                     |
| <b>Dilution Range</b>      | WB 1:500-1:2000                                           |
| <b>Formulation</b>         | PBS containing 0.02% Sodium Azide, 50% Glycerol, pH7.3.   |
| <b>Isotype</b>             | IgG                                                       |
| <b>Storage Instruction</b> | Store in a freezer at -20°C and avoid freeze-thaw cycles. |

### TARGET INFORMATION

|                           |                                                                                                                         |
|---------------------------|-------------------------------------------------------------------------------------------------------------------------|
| <b>Gene ID</b>            | <a href="#">151556</a>                                                                                                  |
| <b>Gene Symbol</b>        | <a href="#">GPR155</a>                                                                                                  |
| <b>Uniprot ID</b>         | <a href="#">GP155_HUMAN</a>                                                                                             |
| <b>Immunogen</b>          | Recombinant fusion protein containing a sequence corresponding to amino acids 721-870 of human GPR155 (NP_001028217.1). |
| <b>Immunogen Region</b>   | 721-870                                                                                                                 |
| <b>Specificity</b>        |                                                                                                                         |
| <b>Immunogen Sequence</b> |                                                                                                                         |



Western blot analysis of extracts of various cell lines, using GPR155 antibody (STJ118613) at 1:1000 dilution. Secondary antibody: HRP Goat Anti-rabbit IgG (H+L) at 1:10000 dilution. Lysates/proteins: 25µg per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL Basic Kit. Exposure time: 5s.

This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.  
St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081