

Anti-SV2B antibody (415-540) (STJ118704)

STJ118704

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-SV2B (415-540) is suitable for use in Western Blot.
Applications	WB
Host/Source	Rabbit
Reactivity	Human, Mouse, Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	WB 1:500-1:2000
Formulation	PBS containing 0.02% Sodium Azide, 50% Glycerol, pH7.3.
Isotype	IgG
Storage Instruction	Store in a freezer at -20°C and avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	9899
Gene Symbol	SV2B
Uniprot ID	SV2B_HUMAN
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 415-540 of human SV2B (NP_055663.1).
Immunogen Region	415-540
Specificity	
Immunogen Sequence	