

Anti-COL17A1 antibody (1300-1400) (STJ118789)

STJ118789

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-COL17A1 (1300-1400) is suitable for use in Western Blot.
Applications	WB
Host/Source	Rabbit
Reactivity	Human, Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	WB 1:500-1:2000
Formulation	PBS containing 0.02% Sodium Azide, 50% Glycerol, pH7.3.
Isotype	IgG
Storage Instruction	Store in a freezer at -20°C and avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	1308
Gene Symbol	COL17A1
Uniprot ID	COHA1_HUMAN
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 1300-1400 of human Collagen XVII/Collagen XVII/COL17A1 (NP_000485.3).
Immunogen Region	1300-1400
Specificity	
Immunogen Sequence	