

Anti-FURIN antibody (C-Term) (STJ118839)

STJ118839

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-FURIN (C-Term) is suitable for use in Western Blot and Immunoprecipitation.
Applications	WB, IP
Host/Source	Rabbit
Reactivity	Human, Mouse

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	WB 1:500-1:2000 IP 1:50-1:200
Formulation	PBS containing 0.02% Sodium Azide, 50% Glycerol, pH7.3.
Isotype	IgG
Storage Instruction	Store in a freezer at -20°C and avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	5045
Gene Symbol	FURIN
Uniprot ID	FURIN_HUMAN
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 700 to the C-terminus of human Furin (NP_002560.1).
Immunogen Region	C-Term
Specificity	
Immunogen Sequence	