

Anti-GFI1B antibody (50-150) (STJ118880)

STJ118880

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-GFI1B (50-150) is suitable for use in Western Blot.
Applications	WB
Host/Source	Rabbit
Reactivity	Mouse

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	WB 1:500-1:2000
Formulation	PBS containing 0.02% Sodium Azide, 50% Glycerol, pH7.3.
Isotype	IgG
Storage Instruction	Store in a freezer at -20°C and avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	8328
Gene Symbol	GFI1B
Uniprot ID	GFI1B_HUMAN
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 50-150 of human GFI1B (NP_001128503.1).
Immunogen Region	50-150
Specificity	
Immunogen Sequence	