

Anti-TM9SF2 antibody (29-300) (STJ118895)

STJ118895

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-TM9SF2 (29-300) is suitable for use in Western Blot.
Applications	WB
Host/Source	Rabbit
Reactivity	Human, Mouse

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	WB 1:500-1:2000
Formulation	PBS containing 0.02% Sodium Azide, 50% Glycerol, pH7.3.
Isotype	IgG
Storage Instruction	Store in a freezer at -20°C and avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	9375
Gene Symbol	TM9SF2
Uniprot ID	TM9SF2_HUMAN
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 29-300 of human TM9SF2 (NP_004791.1).
Immunogen Region	29-300
Specificity	
Immunogen Sequence	