

Anti-SPINK5 antibody (895-1094) (STJ118916)

STJ118916

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-SPINK5 (895-1094) is suitable for use in Western Blot.
Applications	WB
Host/Source	Rabbit
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	WB 1:500-1:2000
Formulation	PBS containing 0.02% Sodium Azide, 50% Glycerol, pH7.3.
Isotype	IgG
Storage Instruction	Store in a freezer at -20°C and avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	11005
Gene Symbol	SPINK5
Uniprot ID	ISK5_HUMAN
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 895-1094 of human SPINK5 (NP_001121170.1).
Immunogen Region	895-1094
Specificity	
Immunogen Sequence	