

## Anti-RHBG antibody (1-100) (STJ118985)

STJ118985

### GENERAL INFORMATION

|                          |   |
|--------------------------|---|
| <b>Product Type</b>      | Primary antibodies  |
| <b>Short Description</b> | Rabbit polyclonal antibody anti-RHBG (1-100) is suitable for use in Western Blot. |
| <b>Applications</b>      | WB  |
| <b>Host/Source</b>       | Rabbit  |
| <b>Reactivity</b>        | Mouse, Rat  |

### PRODUCT PROPERTIES

|                            |   |
|----------------------------|---|
| <b>Clonality</b>           | Polyclonal  |
| <b>Clone ID</b>            |   |
| <b>Concentration</b>       |   |
| <b>Conjugation</b>         | Unconjugated  |
| <b>Purification</b>        | Affinity purification                                     |
| <b>Dilution Range</b>      | WB 1:500-1:2000   |
| <b>Formulation</b>         | PBS containing 0.02% Sodium Azide, 50% Glycerol, pH7.3.   |
| <b>Isotype</b>             | IgG   |
| <b>Storage Instruction</b> | Store in a freezer at -20°C and avoid freeze-thaw cycles. |

### TARGET INFORMATION

|                           |   |
|---------------------------|---|
| <b>Gene ID</b>            | 57127   |
| <b>Gene Symbol</b>        | RHBG  |
| <b>Uniprot ID</b>         | RHBG_HUMAN  |
| <b>Immunogen</b>          | A synthetic peptide corresponding to a sequence within amino acids 1-100 of human RHBG (NP_065140.3). |
| <b>Immunogen Region</b>   | 1-100   |
| <b>Specificity</b>        |   |
| <b>Immunogen Sequence</b> |   |