

Anti-RAB39B antibody (1-213) (STJ119018)

STJ119018

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-RAB39B (1-213) is suitable for use in Western Blot.
Applications	WB
Host/Source	Rabbit
Reactivity	Mouse

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	WB 1:500-1:2000
Formulation	PBS containing 0.02% Sodium Azide, 50% Glycerol, pH7.3.
Isotype	IgG
Storage Instruction	Store in a freezer at -20°C and avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	116442
Gene Symbol	RAB39B
Uniprot ID	RB39B_HUMAN
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-213 of human RAB39B (NP_741995.1).
Immunogen Region	1-213
Specificity	
Immunogen Sequence	