

Anti-TDRD9 antibody (300-600) (STJ119019)

STJ119019

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-TDRD9 (300-600) is suitable for use in Western Blot.
Applications	WB
Host/Source	Rabbit
Reactivity	Mouse, Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	WB 1:500-1:2000
Formulation	PBS containing 0.02% Sodium Azide, 50% Glycerol, pH7.3.
Isotype	IgG
Storage Instruction	Store in a freezer at -20°C and avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	122402
Gene Symbol	TDRD9
Uniprot ID	TDRD9_HUMAN
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 300-600 of human TDRD9 (NP_694591.2).
Immunogen Region	300-600
Specificity	
Immunogen Sequence	