

Anti-RNF152 antibody (100-200) (STJ119036)

STJ119036

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-RNF152 (100-200) is suitable for use in Western Blot.
Applications	WB
Host/Source	Rabbit
Reactivity	Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	WB 1:500-1:2000
Formulation	PBS containing 0.02% Sodium Azide, 50% Glycerol, pH7.3.
Isotype	IgG
Storage Instruction	Store in a freezer at -20°C and avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	220441
Gene Symbol	RNF152
Uniprot ID	RN152_HUMAN
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 100-200 of human RNF152 (NP_775828.1).
Immunogen Region	100-200
Specificity	
Immunogen Sequence	