

Anti-TERT antibody (C-Term) (STJ119052)

STJ119052

GENERAL INFORMATION


Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-TERT (C-Term) is suitable for use in Western Blot and Immunohistochemistry.
Applications	WB, IHC
Host/Source	Rabbit
Reactivity	Human, Mouse

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	WB 1:500-1:2000 IHC 1:50-1:200
Formulation	PBS containing 0.02% Sodium Azide, 50% Glycerol, pH7.3.
Isotype	IgG
Storage Instruction	Store in a freezer at -20°C and avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	7015
Gene Symbol	TERT
Uniprot ID	TERT_HUMAN
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 1000 to the C-terminus of human TERT (NP_937983.2).
Immunogen Region	C-Term
Specificity	
Immunogen Sequence	



Immunohistochemistry of paraffin-embedded human breast cancer using TERT antibody (STJ119052) at dilution of 1:100 (40x lens).