

Anti-ESR1 antibody (1-190) (STJ119086)

STJ119086

GENERAL INFORMATION

| | |
|--------------------------|--|
| Product Type | Primary antibodies |
| Short Description | Rabbit polyclonal antibody anti-ESR1 (1-190) is suitable for use in Western Blot and Immunofluorescence. |
| Applications | WB, IF |
| Host/Source | Rabbit |
| Reactivity | Human, Mouse, Rat |

PRODUCT PROPERTIES

| | |
|-----------------------|---|
| Clonality | Polyclonal |
| Clone ID | |
| Concentration | |
| Conjugation | Unconjugated |
| Purification | Affinity purification |
| Dilution Range | WB 1:1000-1:4000 IF 1:50-1:200 |
| Formulation | PBS containing 0.02% Sodium Azide, 50% Glycerol, pH7.3. |
| Isotype | IgG |
| Storage | Store in a freezer at -20°C and avoid freeze-thaw cycles. |
| Instruction | |

TARGET INFORMATION

| | |
|---------------------------|---|
| Gene ID | 2099 |
| Gene Symbol | ESR1 |
| Uniprot ID | ESR1_HUMAN |
| Immunogen | Recombinant fusion protein containing a sequence corresponding to amino acids 1-190 of human Estrogen Receptor alpha (ESR Alpha) (NP_000116.2). |
| Immunogen Region | 1-190 |
| Specificity | |
| Immunogen Sequence | |