

Anti-PLA2G5 antibody (1-138) (STJ119123)

STJ119123

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-PLA2G5 (1-138) is suitable for use in Immunofluorescence.
Applications	IF
Host/Source	Rabbit
Reactivity	Mouse

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	IF 1:50-1:200
Formulation	PBS containing 0.02% Sodium Azide, 50% Glycerol, pH7.3.
Isotype	IgG
Storage Instruction	Store in a freezer at -20°C and avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	5322
Gene Symbol	PLA2G5
Uniprot ID	PA2G5_HUMAN
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-138 of human PLA2G5 (NP_000920.1).
Immunogen Region	1-138
Specificity	
Immunogen Sequence	