

Anti-CD70 antibody (STJ119187)

STJ119187

GENERAL INFORMATION


Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-CD70 is suitable for use in Immunofluorescence.
Applications	IF
Host/Source	Rabbit
Reactivity	Mouse, Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	IF 1:50-1:200
Formulation	PBS containing 0.02% Sodium Azide, 50% Glycerol, pH7.3.
Isotype	IgG
Storage Instruction	Store in a freezer at -20°C and avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	970
Gene Symbol	CD70
Uniprot ID	CD70_HUMAN
Immunogen	A synthetic peptide of human CD70.
Immunogen Region	
Specificity	
Immunogen Sequence	



Immunofluorescence analysis of rat spleen using CD70 antibody (STJ119187) at dilution of 1:100. Blue: DAPI for nuclear staining.