

Anti-TAX1BP1 antibody (709-789) (STJ119281)

STJ119281

GENERAL INFORMATION


Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-TAX1BP1 (709-789) is suitable for use in Immunohistochemistry.
Applications	IHC
Host/Source	Rabbit
Reactivity	Human, Rat

PRODUCT PROPERTIES


Clonality	Polyclonal
Clone ID	
Concentration	
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	IHC 1:50-1:200
Formulation	PBS containing 0.02% Sodium Azide, 50% Glycerol, pH7.3.
Isotype	IgG
Storage	Store in a freezer at -20°C and avoid freeze-thaw cycles.
Instruction	

TARGET INFORMATION

Gene ID	8887
Gene Symbol	TAX1BP1
Uniprot ID	TAXB1_HUMAN
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 709-789 of human TRAF6BP/TRAF6BP/TAX1BP1 (NP_006015.4).
Immunogen Region	709-789
Specificity	
Immunogen Sequence	



Immunohistochemistry of paraffin-embedded human lung cancer using TRAF6BP/TRAF6BP/TAX1BP1 antibody (STJ119281) at dilution of 1:100 (40x lens).



Immunohistochemistry of paraffin-embedded human breast cancer using TRAF6BP/TRAF6BP/TAX1BP1 antibody (STJ119281) at dilution of 1:100 (40x lens).

This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.
St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081