

Anti-ASIC4 antibody (350-460) (STJ119332)

STJ119332

GENERAL INFORMATION


Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-ASIC4 (350-460) is suitable for use in Immunohistochemistry.
Applications	IHC
Host/Source	Rabbit
Reactivity	Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	IHC 1:50-1:200
Formulation	PBS containing 0.02% Sodium Azide, 50% Glycerol, pH7.3.
Isotype	IgG
Storage Instruction	Store in a freezer at -20°C and avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	55515
Gene Symbol	ASIC4
Uniprot ID	ASIC4_HUMAN
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 350-460 of human ASIC4 (NP_878267.2).
Immunogen Region	350-460
Specificity	
Immunogen Sequence	



Immunohistochemistry of paraffin-embedded rat brain using ASIC4 antibody (STJ119332) at dilution of 1:100 (40x lens).

This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.
St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081