

Anti-RNF14 antibody (125-474) (STJ119441)

STJ119441

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-RNF14 (125-474) is suitable for use in Western Blot.
Applications	WB
Host/Source	Rabbit
Reactivity	Human, Mouse, Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	WB 1:500-1:2000
Formulation	PBS containing 0.02% Sodium Azide, 50% Glycerol, pH7.3.
Isotype	IgG
Storage Instruction	Store in a freezer at -20°C and avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	9604
Gene Symbol	RNF14
Uniprot ID	RNF14_HUMAN
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 125-474 of human RNF14 (NP_004281.1).
Immunogen Region	125-474
Specificity	
Immunogen Sequence	