

Anti-DHRS7 antibody (STJ119469)

STJ119469

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-DHRS7 is suitable for use in Western Blot.
Applications	WB
Host/Source	Rabbit
Reactivity	Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	WB 1:500-1:1000
Formulation	PBS containing 0.02% Sodium Azide, pH7.3.
Isotype	IgG
Storage Instruction	Store in a freezer at -20°C and avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	51635
Gene Symbol	DHRS7
Uniprot ID	DHRS7_HUMAN
Immunogen	A synthetic peptide of human DHRS7
Immunogen Region	
Specificity	
Immunogen Sequence	