

## Anti-TAOK3 antibody (STJ119471)

STJ119471

### GENERAL INFORMATION

|                          |  |
|--------------------------|--|
| <b>Product Type</b>      | Primary antibodies   |
| <b>Short Description</b> | Rabbit polyclonal antibody anti-TAOK3 is suitable for use in Western Blot. |
| <b>Applications</b>      | WB   |
| <b>Host/Source</b>       | Rabbit   |
| <b>Reactivity</b>        | Mouse  |

### PRODUCT PROPERTIES

|                            |   |
|----------------------------|---|
| <b>Clonality</b>           | Polyclonal  |
| <b>Clone ID</b>            |   |
| <b>Concentration</b>       |   |
| <b>Conjugation</b>         | Unconjugated  |
| <b>Purification</b>        | Affinity purification                                     |
| <b>Dilution Range</b>      | WB 1:500-1:2000   |
| <b>Formulation</b>         | PBS containing 0.02% Sodium Azide, 50% Glycerol, pH7.3.   |
| <b>Isotype</b>             | IgG   |
| <b>Storage Instruction</b> | Store in a freezer at -20°C and avoid freeze-thaw cycles. |

### TARGET INFORMATION

|                           |                                     |
|---------------------------|-------------------------------------|
| <b>Gene ID</b>            | 51347                               |
| <b>Gene Symbol</b>        | TAOK3                               |
| <b>Uniprot ID</b>         | TAOK3_HUMAN                         |
| <b>Immunogen</b>          | Recombinant protein of human TAOK3. |
| <b>Immunogen Region</b>   |                                     |
| <b>Specificity</b>        |                                     |
| <b>Immunogen Sequence</b> |                                     |