

Anti-PDGFA antibody (87-211) (STJ119551)

STJ119551

GENERAL INFORMATION


Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-PDGFA (87-211) is suitable for use in Immunohistochemistry.
Applications	IHC
Host/Source	Rabbit
Reactivity	Human, Mouse

PRODUCT PROPERTIES


Clonality	Polyclonal
Clone ID	
Concentration	
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	IHC 1:50-1:200
Formulation	PBS containing 0.02% Sodium Azide, 50% Glycerol, pH7.3.
Isotype	IgG
Storage Instruction	Store in a freezer at -20°C and avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	5154
Gene Symbol	PDGFA
Uniprot ID	PDGFA_HUMAN
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 87-211 of human PDGFA (NP_002598.4).
Immunogen Region	87-211
Specificity	
Immunogen Sequence	



Immunohistochemistry of paraffin-embedded mouse kidney using PDGFA antibody (STJ119551) at dilution of 1:100 (40x lens).



Immunohistochemistry of paraffin-embedded Human esophagus using PDGFA antibody (STJ119551) at dilution of 1:100 (40x lens).

This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.
St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081