

Anti-DPY30 antibody (1-99) (STJ119825)

STJ119825

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-DPY30 (1-99) is suitable for use in Immunohistochemistry.
Applications	IHC
Host/Source	Rabbit
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	IHC 1:50-1:100
Formulation	PBS containing 0.02% Sodium Azide, 50% Glycerol, pH7.3.
Isotype	IgG
Storage Instruction	Store in a freezer at -20°C and avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	84661
Gene Symbol	DPY30
Uniprot ID	DPY30_HUMAN
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-99 of human DPY30 (NP_115963.1).
Immunogen Region	1-99
Specificity	
Immunogen Sequence	