

Anti-KSR2 antibody (780-890) (STJ119878)

STJ119878

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-KSR2 (780-890) is suitable for use in Western Blot.
Applications	WB
Host/Source	Rabbit
Reactivity	Human, Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	WB 1:500-1:2000
Formulation	PBS containing 0.02% Sodium Azide, 50% Glycerol, pH7.3.
Isotype	IgG
Storage Instruction	Store in a freezer at -20°C and avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	283455
Gene Symbol	KSR2
Uniprot ID	KSR2_HUMAN
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 780-890 of human KSR2 (NP_775869.3).
Immunogen Region	780-890
Specificity	
Immunogen Sequence	