

Anti-FoxO1/FoxO3a/FoxO4 antibody (STJ119956)

STJ119956

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-FoxO1/FoxO3a/FoxO4 is suitable for use in Western Blot and Immunohistochemistry.
Applications	WB, IHC
Host/Source	Rabbit
Reactivity	Human, Mouse, Rat


PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	WB 1:500-1:2000 IHC 1:50-1:200
Formulation	PBS containing 0.02% Sodium Azide, 50% Glycerol, pH7.3.
Isotype	IgG
Storage Instruction	Store in a freezer at -20°C and avoid freeze-thaw cycles.


TARGET INFORMATION

Gene ID
[2308](#)
[4303](#)
[2309](#)
[FOXO1](#)
[FOXO4](#)
[FOXO1_HUMAN](#)
[FOXO4_HUMAN](#)
[FOXO3_HUMAN](#)


Immunogen Recombinant protein of human FOXO1/FOXO3A/FOXO4.
Immunogen Region
Immunogen Specificity
Immunogen Sequence




Western blot analysis of extracts of various cell lines, using FOXO1/FOXO3A/FOXO4 antibody (STJ119956) at 1:1000 dilution. Secondary antibody: HRP Goat Anti-rabbit IgG (H+L) at 1:10000 dilution. Lysates/proteins: 25µg per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL Basic Kit. Exposure time: 180s.



Immunohistochemistry of paraffin-embedded rat brain using FOXO1/FOXO3A/FOXO4 rabbit polyclonal antibody (STJ119956) at dilution of 1:100 (40x lens). Perform high pressure antigen retrieval with 10 mM citrate buffer pH 6.0 before commencing with immunohistochemistry staining protocol.



Immunohistochemistry of paraffin-embedded human colon carcinoma using FOXO1/FOXO3A/FOXO4 rabbit polyclonal antibody (STJ119956) at dilution of 1:100 (40x lens). Perform high pressure antigen retrieval with 10 mM citrate buffer pH 6.0 before commencing with immunohistochemistry staining protocol.



Immunohistochemistry of paraffin-embedded mouse spinal cord using FOXO1/FOXO3A/FOXO4 rabbit polyclonal antibody (STJ119956) at dilution of 1:100 (40x lens). Perform high pressure antigen retrieval with 10 mM citrate buffer pH 6.0 before commencing with immunohistochemistry staining protocol.

This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.
 St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081