

Anti-Pro--BDNF antibody (STJ13100048)

STJ13100048

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Nz White Rabbit polyclonal antibody anti-Pro--BDNF is suitable for use in Immunohistochemistry and Western Blot research applications.
Applications	IHC, WB
Host/Source	NZ White Rabbit
Reactivity	Mouse, Rat, Human

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	
Conjugation	Unconjugated
Purification	IgG
Dilution Range	A concentration 10-50 µg/ml is recommended. The optimal concentration should be determined by the end user.
Formulation	Shipped as lyophilised. Reconstitute in 500 µl of sterile water. Centrifuge to remove any insoluble material.
Isotype	IgG
Storage	Maintain the lyophilised/reconstituted antibodies frozen at -20°C for long term storage and refrigerated at 2-8°C for a shorter term.
Instruction	When reconstituting, glycerol (1:1) may be added for an additional stability. Avoid freeze and thaw cycles.

TARGET INFORMATION

Gene ID	627
Gene Symbol	BDNF
Uniprot ID	BDNF_HUMAN
Immunogen	A synthetic peptide from mouse pro BDNF conjugated to blue carrier protein was used as the antigen.
Immunogen Region	
Specificity	Specific for proBDNF.
Immunogen Sequence	