

Anti-BDNF antibody (STJ13100055)

STJ13100055

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Sheep polyclonal antibody anti-BDNF is suitable for use in Immunohistochemistry and Western Blot research applications.
Applications	IHC, WB
Host/Source	Sheep
Reactivity	Human, Rat, Mouse, Marmoset

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	
Conjugation	Unconjugated
Purification	Affinity purified
Dilution Range	A concentration of 1-5 µg/ml is recommended. The optimal concentration should be determined by the end user. Not yet tested in other applications.
Formulation	Shipped as lyophilised. Reconstitute in 100 µl of sterile water. Centrifuge to remove any insoluble material.
Isotype	
Storage	Maintain the lyophilised/reconstituted antibodies frozen at -20°C for long term storage and refrigerated at 2-8°C for a shorter term.
Instruction	When reconstituting, glycerol (1:1) may be added for an additional stability. Avoid freeze and thaw cycles.

TARGET INFORMATION

Gene ID	627
Gene Symbol	BDNF
Uniprot ID	BDNF_HUMAN
Immunogen	A synthetic peptide from human mature BDNF conjugated to blue carrier protein was used as the antigen. The peptide is homologous in many other species including rat, mouse, zebra fish and xenopus.
Immunogen Region	
Specificity	Specific for mature BDNF.
Immunogen Sequence	