

Anti-BDKRB2 antibody (C-Term) (STJ13100059)

STJ13100059

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Nz White Rabbit polyclonal antibody anti-BDKRB2 (C-Term) is suitable for use in Immunohistochemistry and Western Blot research applications.
Applications	IHC, WB
Host/Source	NZ White Rabbit
Reactivity	Mouse, Rat, Marmoset

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	
Conjugation	Unconjugated
Purification	IgG
Dilution Range	A concentration of 10-50 µg/ml is recommended. The optimal concentration should be determined by the end user. Not yet tested in other applications.
Formulation	Shipped as lyophilised. Reconstitute in 500 µl of sterile water. Centrifuge to remove any insoluble material.
Isotype	IgG
Storage	Maintain the lyophilised/reconstituted antibodies frozen at -20°C for long term storage and refrigerated at 2-8°C for a shorter term.
Instruction	When reconstituting, glycerol (1:1) may be added for an additional stability. Avoid freeze and thaw cycles.

TARGET INFORMATION

Gene ID	12062
Gene Symbol	Bdkrb2
Uniprot ID	BKRB2_MOUSE
Immunogen	A synthetic peptide from the c-terminal region of mouse BDKRB2 conjugated to an immunogenic carrier protein was used as the antigen.
Immunogen Region	C-Term
Specificity	Specific for BDKRB2.
Immunogen Sequence	