

Anti-C1QTNF5 antibody (Internal) (STJ13100082)

STJ13100082

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Nz White Rabbit polyclonal antibody anti-CTRP5 (Internal) is suitable for use in Immunohistochemistry and Western Blot research applications.
Applications	IHC, WB
Host/Source	NZ White Rabbit
Reactivity	Human, Rat, Mouse

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	
Conjugation	Unconjugated
Purification	Whole serum
Dilution Range	A dilution of 1:300 to 1:2000 is recommended. The optimal dilution should be determined by the end user. Not yet tested in other applications.
Formulation	Shipped as lyophilised. Reconstitute in 100 µl of sterile water. Centrifuge to remove any insoluble material.
Isotype	IgG
Storage	Maintain the lyophilised/reconstituted antibodies frozen at -20°C for long term storage and refrigerated at 2-8°C for a shorter term.
Instruction	When reconstituting, glycerol (1:1) may be added for an additional stability. Avoid freeze and thaw cycles.

TARGET INFORMATION

Gene ID	114902
Gene Symbol	C1QTNF5
Uniprot ID	C1QT5_HUMAN
Immunogen	A synthetic peptide from internal region of human CTRP5 conjugated to an immunogenic carrier protein was used as the antigen. The antigen is homologous in rat and mouse.
Immunogen Region	Internal
Specificity	Specific for CTRP5.
Immunogen Sequence	