

Anti-OR56A4 antibody (STJ13100247)

STJ13100247

GENERAL INFORMATION

| | |
|--------------------------|---|
| Product Type | Primary antibodies |
| Short Description | Nz White Rabbit polyclonal antibody anti-OR56A4 is suitable for use in Immunohistochemistry and Western Blot research applications. |
| Applications | IHC, WB |
| Host/Source | NZ White Rabbit |
| Reactivity | Human |

PRODUCT PROPERTIES

| | |
|-----------------------|---|
| Clonality | Polyclonal |
| Clone ID | |
| Concentration | |
| Conjugation | Unconjugated |
| Purification | Whole serum |
| Dilution Range | A dilution of 1:300 to 1:2000 is recommended. The optimal dilution should be determined by the end user. |
| Formulation | Shipped as lyophilised. Reconstitute in 100 µl of sterile water. Centrifuge to remove any insoluble material. |
| Isotype | IgG |
| Storage | Maintain the lyophilised/reconstituted antibodies frozen at -20°C for long term storage and refrigerated at 2-8°C for a shorter term. |
| Instruction | When reconstituting, glycerol (1:1) may be added for an additional stability. Avoid freeze and thaw cycles. |

TARGET INFORMATION

| | |
|---------------------------|---|
| Gene ID | 120793 |
| Gene Symbol | OR56A4 |
| Uniprot ID | O56A4_HUMAN |
| Immunogen | A synthetic peptide from human olfactory receptor family (OR56A4) conjugated to blue carrier protein was used as the antigen. |
| Region | |
| Specificity | Specific for OR56A4 |
| Immunogen Sequence | |