

Anti-RAB6B antibody (C-Term) (STJ13100309)

STJ13100309

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Nz White Rabbit polyclonal antibody anti-RAB6B (C-Term) is suitable for use in Immunohistochemistry and Western Blot research applications.
Applications	IHC, WB
Host/Source	NZ White Rabbit
Reactivity	Human, Mouse

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	
Conjugation	Unconjugated
Purification	Whole serum
Dilution Range	Use at a dilution of 1:300 to 1:2000. The optimal dilution should be determined by the end user. Not yet tested in other applications.
Formulation	Shipped as lyophilised. Reconstitute in 100 µl of sterile water. Centrifuge to remove any insoluble material.
Isotype	IgG
Storage	Maintain the lyophilised/reconstituted antibodies frozen at -20°C for long term storage and refrigerated at 2-8°C for a shorter term.
Instruction	When reconstituting, glycerol (1:1) may be added for an additional stability. Avoid freeze and thaw cycles.

TARGET INFORMATION

Gene ID	51560
Gene Symbol	RAB6B
Uniprot ID	RAB6B_HUMAN
Immunogen	A synthetic peptide from c-terminal region of human Ras-related protein Rab-6B (RAB6B) conjugated to an immunogenic carrier protein has been used as the antigen.
Immunogen Region	C-Term
Specificity	Specific for RAB6B.
Immunogen Sequence	