

Anti-TRPA1 antibody (STJ13100386)

STJ13100386

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Nz White Rabbit polyclonal antibody anti-TRPA1 is suitable for use in Immunohistochemistry and Western Blot research applications.
Applications	IHC, WB
Host/Source	NZ White Rabbit
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	
Conjugation	Unconjugated
Purification	IgG
Dilution Range	Use at a concentration of 10-50 µg/ml. The optimal working concentration should be determined by the end user.
Formulation	Shipped as lyophilised. Reconstitute in 500 µl of sterile water. Centrifuge to remove any insoluble material.
Isotype	IgG
Storage	Maintain the lyophilised/reconstituted antibodies frozen at -20°C for long term storage and refrigerated at 2-8°C for a shorter term.
Instruction	When reconstituting, glycerol (1:1) may be added for an additional stability. Avoid freeze and thaw cycles.

TARGET INFORMATION

Gene ID	8989
Gene Symbol	TRPA1
Uniprot ID	TRPA1_HUMAN
Immunogen	A synthetic peptide from human TRPA1 conjugated to blue carrier protein was used as the antigen.
Immunogen Region	
Specificity	Specific for TRPA1.
Immunogen Sequence	