

## Anti-TRPV4 antibody (300-400) (STJ13100422)

STJ13100422

### GENERAL INFORMATION

<b>Product Type</b>	Primary antibodies
<b>Short Description</b>	Nz White Rabbit polyclonal antibody anti-TRPV4 (300-400) is suitable for use in Immunohistochemistry and Western Blot research applications.
<b>Applications</b>	IHC, WB
<b>Host/Source</b>	NZ White Rabbit
<b>Reactivity</b>	Mouse, Human, Rat, Marmoset

### PRODUCT PROPERTIES

<b>Clonality</b>	Polyclonal
<b>Clone ID</b>	
<b>Concentration</b>	
<b>Conjugation</b>	Unconjugated
<b>Purification</b>	IgG
<b>Dilution Range</b>	A concentration of 10-50 µg/ml is recommended. The optimal concentration should be determined by the end user. Not yet tested in other application.
<b>Formulation</b>	Shipped as lyophilised. Reconstitute in 500 µl of sterile water. Centrifuge to remove any insoluble material.
<b>Isotype</b>	IgG
<b>Storage</b>	Maintain the lyophilised/reconstituted antibodies frozen at -20°C for long term storage and refrigerated at 2-8°C for a shorter term.
<b>Instruction</b>	When reconstituting, glycerol (1:1) may be added for an additional stability. Avoid freeze and thaw cycles.

### TARGET INFORMATION

<b>Gene ID</b>	<a href="#">59341</a>
<b>Gene Symbol</b>	<a href="#">TRPV4</a>
<b>Uniprot ID</b>	<a href="#">TRPV4_HUMAN</a>
<b>Immunogen</b>	A synthetic peptide from aa region 300-400 of mouse TRPV4 conjugated to an immunogenic carrier protein was used as the antigen. The peptide is homologous in rat and human.
<b>Immunogen Region</b>	300-400
<b>Specificity</b>	Specific for TRPV4.
<b>Immunogen Sequence</b>	