

Anti-VAMP1, 2, 3 antibody (25-65) (STJ13100474)

STJ13100474

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Nz White Rabbit polyclonal antibody anti-VAMP1, 2, 3 (25-65) is suitable for use in Immunohistochemistry and Western Blot research applications.
Applications	IHC, WB
Host/Source	NZ White Rabbit
Reactivity	Human, Rat, Mouse

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	
Conjugation	Unconjugated
Purification	Whole serum
Dilution Range	A working concentration of 10-50 µg/ml is recommended. The optimal concentration should be determined by the end user.
Formulation	Shipped as lyophilised. Reconstitute in 100 µl of sterile water. Centrifuge to remove any insoluble material.
Isotype	IgG
Storage	Maintain the lyophilised/reconstituted antibodies frozen at -20°C for long term storage and refrigerated at 2-8°C for a shorter term.
Instruction	When reconstituting, glycerol (1:1) may be added for an additional stability. Avoid freeze and thaw cycles.

TARGET INFORMATION

Gene ID	6843 6844 9341 VAMP1 VAMP2 VAMP1_HUMAN VAMP2_HUMAN VAMP3_HUMAN
Immunogen	A synthetic peptide from aa region 25-65 of human VAMP 3 conjugated to an immunogenic carrier protein was used as the antigen.
Immunogen Region	25-65
Specificity	Specific for VAMP1, 2, 3.
Immunogen Sequence	