

Anti-TGN38 antibody (STJ13100507)

STJ13100507

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Guinea Pig polyclonal antibody anti-TGN38 is suitable for use in Immunohistochemistry and Western Blot research applications.
Applications	IHC, WB
Host/Source	Guinea Pig
Reactivity	Rat, Mouse, Human, Marmoset

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	
Conjugation	Unconjugated
Purification	Whole serum
Dilution Range	A dilution of 1: 500 is recommended. The optimal dilution should be determined by the end user. Not yet tested in other applications.
Formulation	Shipped as lyophilised. Reconstitute in 100 ul of sterile water. Centrifuge to remove any insoluble material.
Isotype	
Storage	Maintain the lyophilised/reconstituted antibodies frozen at -20°C for long term storage and refrigerated at 2-8°C for a shorter term.
Instruction	When reconstituting, glycerol (1:1) may be added for an additional stability. Avoid freeze and thaw cycles.

TARGET INFORMATION

Gene ID	22135
Gene Symbol	Tgoln2
Uniprot ID	TGON2_MOUSE
Immunogen	A synthetic peptide from mouse TGN38 conjugated to blue carrier protein was used as the antigen.
Immunogen Region	
Specificity	Specific for TGN38 (or TGN46, TGN48, TGN51 in human).
Immunogen Sequence	