

Anti-SARS-CoV-2 NSP5 antibody (1-180) (STJ140186)

STJ140186

GENERAL INFORMATION


Product Type	Primary antibodies
Short Description	Goat polyclonal antibody anti-SARS-CoV-2 NSP5 (1-180) is suitable for use in Western Blot research applications.
Applications	WB
Host/Source	Goat
Reactivity	SARS-CoV-2

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	3 mg/mL
Conjugation	Unconjugated
Purification	Epitope affinity purified
Dilution Range	WB:1:500-1:2000
Formulation	PBS, 20% glycerol and 0.05% sodium azide
Isotype	IgG
Storage	For continuous use, store at 2-8 C for one-two days. For extended storage, store in -20 C freezer. Working dilution samples should be discarded if not used within 12 hours.
Instruction	

TARGET INFORMATION

Gene ID	
Gene Symbol	
Uniprot ID	
Immunogen	Affinity purified recombinant fusion protein using the N-terminal of Nsp5 (residues 1 to 180) and produced in E. coli.
Immunogen Region	1-180
Specificity	In lysates of transfected cells with the plasmid containing the sequence used, detects the fusion protein by Western blot.
Immunogen Sequence	



Anti-NSP5 antibody at 1/2, 500 dilution; lane with 30 ng of recombinant fusion protein; rabbit polyclonal to goat IgG (HRP) at 1:10000 dilution

This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.
 St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081