

Anti-SARS-CoV-2 ORF3a antibody (160-stop) (STJ140201)

STJ140201

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Goat polyclonal antibody anti-SARS-CoV-2 ORF3a (160-stop) is suitable for use in Western Blot research applications.
Applications	WB
Host/Source	Goat
Reactivity	SARS-CoV-2

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	3 mg/mL
Conjugation	Unconjugated
Purification	Epitope affinity purified
Dilution Range	WB:1:500-1:2000
Formulation	PBS, 20% glycerol and 0.05% sodium azide
Isotype	IgG
Storage	For continuous use, store at 2-8 C for one-two days. For extended storage, store in -20 C freezer. Working dilution samples should be discarded if not used within 12 hours.
Instruction	

TARGET INFORMATION

Gene ID	
Gene Symbol	
Uniprot ID	
Immunogen	Affinity purified fusion recombinant protein using the C-terminal of ORF3a (residues 160 to stop) and produced in E. coli.
Immunogen Region	160-stop
Specificity	In lysates of transfected cells with the plasmid containing the sequence used, detects the fusion protein by Western blot.
Immunogen Sequence	