

Anti-GFP antibody {DyLight®550} (STJ140276)

STJ140276

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Goat polyclonal antibody anti-GFP is suitable for use in Western Blot, Immunofluorescence and Immunohistochemistry research applications.
Applications	WB, IF, IHC-F
Host/Source	Goat
Reactivity	GFP

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	2.5 mg/mL
Conjugation	DyLight®550
Purification	Epitope affinity purified
Dilution Range	WB:1:500-1:5000 IF:1:50-1:1000 IHC-F:1:50-1:1000
Formulation	PBS, 20% glycerol and 0.05% sodium azide
Isotype	IgG
Storage Instruction	For continuous use, store at 2-8 C for one-two days. For extended storage, store in -20 C freezer. Working dilution samples should be discarded if not used within 12 hours.

TARGET INFORMATION

Gene ID	
Gene Symbol	
Uniprot ID	
Immunogen	Purified recombinant peptide produced in E. coli.
Immunogen Region	
Specificity	In 293HEK cells transfected with cds plasmid detects a band of 27 kDa by Western blot. This antibody does not recognize mCherry fluorescent protein.
Immunogen Sequence	