

Mouse beta-42 (Amyloid Beta 1-42) Sandwich ELISA Kit (STJ150004)

STJ150004

GENERAL INFORMATION


| | |
|--------------------------|---|
| Product Type | Kits |
| Description Short | This Mouse beta-42 (Amyloid Beta 1-42) Sandwich ELISA Kit is an in-vitro enzyme-linked immunosorbent assay for the measurement of samples in cell culture supernatant, serum and plasma (EDTA, citrate, heparin). |
| Applications | ELISA |
| Reactivity | Mouse |

PRODUCT PROPERTIES

| | |
|------------------------------|------------------------------|
| Detection Sensitivity | 0.313--20ng/mL |
| Detection Limit | 0.313--20ng/mL |
| Assay Time | |
| Storage | Store at 4°C for 6 months. |
| Instruction | |
| Note | FOR RESEARCH USE ONLY (RUO). |

TARGET INFORMATION

| | |
|---------------------|----------------------------------|
| Sample Types | Serum, Plasma, Biological Fluids |
| Gene | Gh1 |
| Symbol | |
| Uniprot ID | SOMA_MOUSE |



This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.
St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081