

## Mouse C-P (C-Peptide) Sandwich ELISA Kit (STJ150015)

STJ150015

### GENERAL INFORMATION


|                          |   |
|--------------------------|---|
| <b>Product Type</b>      | Kits  |
| <b>Short Description</b> | This Mouse C-P (C-Peptide) Sandwich ELISA Kit is an in-vitro enzyme-linked immunosorbent assay for the measurement of samples in cell culture supernatant, serum and plasma (EDTA, citrate, heparin). |
| <b>Applications</b>      | ELISA   |
| <b>Reactivity</b>        | Mouse   |

### PRODUCT PROPERTIES

|                              |                              |
|------------------------------|------------------------------|
| <b>Detection Sensitivity</b> | 15.625--1000pg/mL            |
| <b>Detection Limit</b>       | 15.625--1000pg/mL            |
| <b>Assay Time</b>            |                              |
| <b>Storage</b>               | Store at 4°C for 6 months.   |
| <b>Instruction</b>           |                              |
| <b>Note</b>                  | FOR RESEARCH USE ONLY (RUO). |

### TARGET INFORMATION

|                     |                                  |
|---------------------|----------------------------------|
| <b>Sample Types</b> | Serum, Plasma, Biological Fluids |
| <b>Gene</b>         | <a href="#">Calca</a>            |
| <b>Symbol</b>       |                                  |
| <b>Uniprot ID</b>   | <a href="#">CALC_MOUSE</a>       |



This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.  
St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081