

Rat PAF (Platelet Activating Factor) Sandwich ELISA Kit (STJ150082)

STJ150082

GENERAL INFORMATION

Product Type	Kits
Short Description	This Rat PAF (Platelet Activating Factor) Sandwich ELISA Kit is an in-vitro enzyme-linked immunosorbent assay for the measurement of samples in cell culture supernatant, serum and plasma (EDTA, citrate, heparin).
Applications	ELISA
Reactivity	Rat

PRODUCT PROPERTIES

Detection Sensitivity	0.156--10ng/mL
Detection Limit	0.156--10ng/mL
Assay Time	
Storage	Store at 4°C for 6 months.
Instruction	
Note	FOR RESEARCH USE ONLY (RUO).

TARGET INFORMATION

Sample Types	Serum, Plasma, Biological Fluids
Gene	Ccl2
Symbol	
Uniprot ID	CCL2_RAT