

Rat VEGFC Sandwich ELISA Kit (STJ150087)

STJ150087

GENERAL INFORMATION


| | |
|--------------------------|---|
| Product Type | Kits |
| Short Description | This Rat VEGFC Sandwich ELISA Kit is an in-vitro enzyme-linked immunosorbent assay for the measurement of samples in cell culture supernatant, serum and plasma (EDTA, citrate, heparin). |
| Applications | ELISA |
| Reactivity | Rat |

PRODUCT PROPERTIES

| | |
|------------------------|------------------------------|
| Detection | 31.25--2000pg/mL |
| Sensitivity | 2000pg/mL |
| Detection Limit | 31.25--2000pg/mL |
| Assay Time | |
| Storage | Store at 4°C for 6 months. |
| Instruction | |
| Note | FOR RESEARCH USE ONLY (RUO). |

TARGET INFORMATION

| | |
|---------------------|----------------------------------|
| Sample Types | Serum, Plasma, Biological Fluids |
| Gene | VEGFC |
| Symbol | |
| Uniprot ID | VEGFC_HUMAN |



This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.
 St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081