

Human OxLDL (Oxidized Low Density Lipoprotein) Sandwich ELISA Kit (STJ150114)

STJ150114

GENERAL INFORMATION


Product Type	Kits
Short Description	This Human OxLDL (Oxidized Low Density Lipoprotein) Sandwich ELISA Kit is an in-vitro enzyme-linked immunosorbent assay for the measurement of samples in cell culture supernatant, serum and plasma (EDTA, citrate, heparin).
Applications	ELISA
Reactivity	Human

PRODUCT PROPERTIES

Detection Limit	31.25--2000pg/mL
Sensitivity	2000pg/mL
Detection Limit	31.25--2000pg/mL
Assay Time	
Storage	Store at 4°C for 6 months.
Instruction	
Note	FOR RESEARCH USE ONLY (RUO).

TARGET INFORMATION

Sample Types	Serum, Plasma, Biological Fluids
Gene	
Symbol	
Uniprot ID	



This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.
St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081