

Mouse ALPL Sandwich ELISA Kit (STJ150162)

STJ150162

GENERAL INFORMATION

| | |
|--------------------------|--|
| Product Type | Kits |
| Short Description | This Mouse ALPL Sandwich ELISA Kit is an in-vitro enzyme-linked immunosorbent assay for the measurement of samples in cell culture supernatant, serum and plasma (EDTA, citrate, heparin). |
| Applications | ELISA |
| Reactivity | Mouse |

PRODUCT PROPERTIES

| | |
|------------------------|------------------------------|
| Detection | 0.156--10ng/mL |
| Sensitivity | 10ng/mL |
| Detection Limit | 0.156--10ng/mL |
| Assay Time | |
| Storage | Store at 4°C for 6 months. |
| Instruction | |
| Note | FOR RESEARCH USE ONLY (RUO). |

TARGET INFORMATION

| | |
|---------------------|----------------------------------|
| Sample Types | Serum, Plasma, Biological Fluids |
| Gene | ALPL |
| Symbol | |
| Uniprot ID | PPBT_HUMAN |