

Mouse PAF (Platelet Activating Factor) Sandwich ELISA Kit (STJ150167)

STJ150167

GENERAL INFORMATION

Product Type	Kits
Short Description	This Mouse PAF (Platelet Activating Factor) Sandwich ELISA Kit is an in-vitro enzyme-linked immunosorbent assay for the measurement of samples in cell culture supernatant, serum and plasma (EDTA, citrate, heparin).
Applications	ELISA
Reactivity	Mouse

PRODUCT PROPERTIES

Detection	0.781--50ng/mL
Sensitivity	50ng/mL
Detection Limit	0.781--50ng/mL
Assay Time	
Storage	Store at 4°C for 6 months.
Instruction	
Note	FOR RESEARCH USE ONLY (RUO).

TARGET INFORMATION

Sample Types	Serum, Plasma, Biological Fluids
Gene	
Symbol	
Uniprot ID	