

Human Bcl-2 (B-cell Leukemia/Lymphoma 2) Sandwich ELISA Kit (STJ150242)

STJ150242

GENERAL INFORMATION


| | |
|---------------------|--|
| Product Type | Kits |
| Description | This Human Bcl-2 (B-cell Leukemia/Lymphoma 2) Sandwich ELISA Kit is an in-vitro enzyme-linked immunosorbent assay for the measurement of samples in cell culture supernatant, serum and plasma (EDTA, citrate, heparin). |
| Applications | ELISA |
| Reactivity | Human |

PRODUCT PROPERTIES

| | |
|------------------------|------------------------------|
| Detection Limit | 0.156--10ng/mL |
| Sensitivity | 10ng/mL |
| Detection Limit | 0.156--10ng/mL |
| Assay Time | |
| Storage | Store at 4°C for 6 months. |
| Instruction | |
| Note | FOR RESEARCH USE ONLY (RUO). |

TARGET INFORMATION

| | |
|---------------------|----------------------------------|
| Sample Types | Serum, Plasma, Biological Fluids |
| Gene | Bcl2 |
| Symbol | |
| Uniprot ID | BCL2_MOUSE |



This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.
St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081