

Rat ICTP (Cross-linked C-Telopeptide of Type I Collagen) Sandwich ELISA Kit (STJ150352)

STJ150352

GENERAL INFORMATION

Product Type	Kits
Short Description	This Rat ICTP (Cross-linked C-Telopeptide of Type I Collagen) Sandwich ELISA Kit is an in-vitro enzyme-linked immunosorbent assay for the measurement of samples in cell culture supernatant, serum and plasma (EDTA, citrate, heparin).
Applications	ELISA
Reactivity	Rat

PRODUCT PROPERTIES

Detection	0.313--20ng/mL
Sensitivity	20ng/mL
Detection Limit	0.313--20ng/mL
Assay Time	
Storage	Store at 4°C for 6 months.
Instruction	
Note	FOR RESEARCH USE ONLY (RUO).

TARGET INFORMATION

Sample Types	Serum, Plasma, Biological Fluids
Gene	
Symbol	
Uniprot ID	